

Philips HeartStart Review Express Connect Product Information



Download and send defibrillator data where it needs to go

HeartStart Review Express Connect is an easy-to-use *wizard* that guides you through the steps of downloading an ECG from a Philips or Laerdal defibrillator, allowing you to:

- view and print the ECG
- save it to a file
- e-mail it to a central data manager or medical director
- erase patient data from the defibrillator's data card or internal memory.

You choose only the steps you want to perform, and Review Express Connect remembers them for the next use.

Review Express Connect is particularly helpful when you simply want to download a case from a defibrillator and e-mail it to a central data manager or medical director for analysis using the more comprehensive HeartStart Event Review or Event Review Pro data management program.

To ensure your continued satisfaction, we support you with expert technical assistance via telephone or e-mail, a knowledge base on the Event Review website, www.medical.philips.com/goto/eventreview, and product updates available for download over the Internet.



Technical Specifications

HeartStart Review Express Connect

This product is provided on CD. An Internet connection is recommended for software updates.

Ordering Product Number

861311 Opt A01

Hardware

Processor

Pentium® processor, 500 MHz or higher

Display

Minimum: Super VGA 800x600 16-bit color

Recommended: 1024x768 screen resolution with 16-bit color

Memory

256 MB minimum

Hard Disk

300 MB of available disk space during software installation and 100 MB minimum for event storage. Hard disk space requirements will vary depending on usage. Variables affecting disk space requirements are the number of cases being archived and archival of audio data. A 15-minute ECG with no audio is approximately 100 KB; with audio, it can be as large as 4.5 MB.

Drive

CD-ROM drive

Datacard Reader

A datacard reader or infrared reader is required in order to transfer defibrillator data to a PC. Visit our website, www.medical.philips.com/goto/eventreview, for details.

Modem

Modem or Internet connection is required to send cases by e-mail

Software

Operating System

Windows 2000 Professional (with SP4) or Windows XP Pro (with SP1 or later)

Browser

Microsoft Internet Explorer 6.0 and later

E-mail

MAPI-compliant e-mail client such as Microsoft Outlook

Defibrillators Supported

Philips & Laerdal Heartstart HSI Family

Philips & Laerdal FRx

Philips & Laerdal Heartstart FR2-series

Philips ForeRunner

Philips HeartStart MRx

Philips HeartStart XL

Philips HeartStart XLT

Laerdal Heartstart FR

Laerdal Heartstart 4000

Call Philips Today

To learn more about HeartStart Data Management Solutions and Philips Medical Systems, contact your local sales representative or visit www.philips.com/heartstart.

**Philips Medical Systems
is part of Royal Philips
Electronics**

www.philips.com/heartstart
medical@philips.com

Philips Medical Systems
3000 Minuteman Road
Andover, MA 01810-1085
+ 1 800 934-7372

© Koninklijke Philips Electronics
N.V. 2005. All rights reserved.
Reproduction in whole or in part
is prohibited without prior written
consent of the copyright holder.

Philips Medical Systems North
America Corporation reserves
the right to make changes in
specifications or to discontinue any
product at any time without notice
or obligation and will not be liable
for any consequences resulting
from the use of this publication.

Philips is a registered trademark of
Koninklijke Philips Electronics N.V.

Printed in The Netherlands.
4522 962 00091/861 * MAR 2005