



ACCESSORIES

for use with

LIFEPAK CR® Plus

AND LIFEPAK EXPRESS® DEFIBRILLATORS



Genuine accessories from Physio-Control

**PHYSIO
CONTROL**

BATTERY AND ELECTRODE OPTIONS

for LIFEPAK CR Plus/LIFEPAK EXPRESS DEFIBRILLATORS



LIFEPAK CR Plus AED

LIFEPAK EXPRESS AED

Thank you for choosing Physio-Control as your partner in helping you save lives and improve patient care. You already know our reputation for quality and reliability in LIFEPAK defibrillators and AEDs. When it comes to LIFEPAK accessories, we know you will not accept anything less. Inside this catalog you will find accessories designed to meet your needs and reliability standards. All the accessories and disposables in this catalog are designed to work with LIFEPAK CR Plus AED and LIFEPAK EXPRESS AEDs. If you have any questions about accessories or disposables, please contact your Physio-Control representative.



Replacement Kit for CHARGE-PAK™ Battery Charger

Includes 2 sets electrodes and 1 battery charger, replacement instructions, and discharger for safe disposal of used CHARGE-PAK.....11403-00001



Replacement Kit for CHARGE-PAK Battery Charger

Includes 1 set electrodes and 1 battery charger, replacement instructions, and discharger for safe disposal of used CHARGE-PAK.

11403-00002



Infant/Child Reduced Energy Defibrillation Electrode Starter Kit

For use only with LIFEPAK® 500 Biphasic AED with pink connector or LIFEPAK 1000 defibrillator, LIFEPAK EXPRESS AED or LIFEPAK CR Plus AED. For use on children less than 8 years of age or less than 55 lbs.

11101-000017



Replacement Infant/Child Reduced Energy Defibrillation Electrodes

11101-000016

CARRYING CASES

for LIFEPAK CR Plus/LIFEPAK EXPRESS DEFIBRILLATORS



Hard-shell, Water-tight Carrying Case11260-00015



Complete Soft Shell Carrying Case.....21300-004576

Replacement Shoulder Strap for Carrying Case.....21300-004579



Replacement Lid.....21300-004913

MISCELLANEOUS ACCESSORIES

for LIFEPAK CR Plus/LIFEPAK EXPRESS DEFIBRILLATORS



AMBU® Res-Cue Mask® First Responder Kit

Includes 2 sets of disposable vinyl gloves, 1 reusable mouth barrier mask with valve and filter, 1 disposable razor, 1 disposable anti-microbial wipe, and 1 pair of trauma scissors. Can be attached to the AED handle.

11998-000321

AMBU Res-Cue Key® First Responder Kit

Includes 2 sets of disposable vinyl gloves, 1 disposable CPR barrier, 1 disposable razor, 1 disposable anti-microbial wipe, and 1 pair of trauma scissors. Can be stored in the AED carrying case **11998-000320**



AED Wall Cabinet with Alarm

Surface mount, rolled edges. Works with LIFEPAK 500, LIFEPAK 1000, LIFEPAK CR Plus, or LIFEPAK EXPRESS defibrillators. Steel finish wall cabinet with white trim. Surface mounted trim style with 7" return.

11220-000079



AED Wall Cabinet with Alarm

Recessed mount, rolled edges. Works with LIFEPAK 500, LIFEPAK 1000, LIFEPAK CR Plus, or LIFEPAK EXPRESS defibrillators. Steel finish wall cabinet with white trim. Recessed mounted trim style with 1.5" return.

11998-000293

AED Wall Cabinet with Alarm

Semi-recessed, rolled edges. Works with LIFEPAK 500, LIFEPAK 1000, LIFEPAK CR Plus, or LIFEPAK EXPRESS defibrillators. Steel finish wall cabinet with white trim. Surface mounted trim style with 3" return.

11998-000292

AED Wall Cabinet with Alarm, Fire Rated

Semi-recessed, rolled edges. Works with LIFEPAK 500, LIFEPAK 1000, LIFEPAK CR Plus, or LIFEPAK EXPRESS defibrillators. Steel finish wall cabinet with white trim with 3" return.

11210-000026

AED Wall Cabinet with Alarm, Fire Rated

Recessed, square edges. Works with LIFEPAK 500, LIFEPAK 1000, LIFEPAK CR Plus, or LIFEPAK EXPRESS defibrillators. Steel finish wall cabinet with white trim. Recessed mounted trim style with 1.5" return.

11210-000027

AED Wall Cabinet with Alarm

Surface mount, rolled edges. Works with LIFEPAK 500, LIFEPAK 1000, LIFEPAK CR Plus, or LIFEPAK EXPRESS defibrillators. Steel, with stainless steel trim. Surface mounted trim style with 7" return.

11220-000076



Wall Mount Bracket for LIFEPAK CR Plus or LIFEPAK EXPRESS AEDs.....

11210-000021

MISCELLANEOUS ACCESSORIES

for LIFEPAK CR Plus/LIFEPAK EXPRESS DEFIBRILLATORS



AED Wall Cabinet with Alarm and Strobe

Surface mount, rolled edges. Works with LIFEPAK 500, LIFEPAK 1000, LIFEPAK CR Plus, or LIFEPAK EXPRESS defibrillators. Steel, with stainless steel trim. Surface mounted trim style with 7" return.

11220-000084

AED Wall Cabinet with Alarm

Semi-recessed, rolled edges. Works with LIFEPAK 500, LIFEPAK 1000, LIFEPAK CR Plus, or LIFEPAK EXPRESS defibrillators. Steel wall cabinet with stainless steel trim. Semi-recessed mounted trim style with 3" return.

11220-000077

AED Wall Cabinet with Alarm

Recessed, rolled edges. Works with LIFEPAK 500, LIFEPAK 1000, LIFEPAK CR Plus, or LIFEPAK EXPRESS defibrillators. Steel finish wall cabinet with stainless steel trim. Recessed mounted trim style with 1.5" return.

11220-000078

AED Wall Cabinet with Alarm and Strobe

Surface mount, rolled edges. Works with LIFEPAK 500, LIFEPAK 1000, LIFEPAK CR Plus, or LIFEPAK EXPRESS defibrillators. Steel finish wall cabinet with white trim. Surface mounted trim style with 7" return.

11220-000083



AED Cabinet Window Replacement Kit

21300-006797

Central Alarm Switch

For use with LIFEPAK CR Plus AED only.

21300-006587



AED Floor Stand Cabinet with Alarm

Works with LIFEPAK 500, LIFEPAK 1000, LIFEPAK CR Plus, or LIFEPAK EXPRESS defibrillators. White epoxy painted steel 52" tall, 14" x 17-3/8" high x 7" deep 11210-000028

AED Floor Stand Cabinet with Alarm

Works with LIFEPAK 500, LIFEPAK 1000, LIFEPAK CR Plus, or LIFEPAK EXPRESS defibrillators. Gray epoxy painted steel. 52" tall, 14" x 17-3/8" high x 7" deep 11210-000029

MISCELLANEOUS ACCESSORIES/TRAINING TOOLS

for LIFEPAK CR Plus/LIFEPAK EXPRESS DEFIBRILLATORS



AED Location Sign

Flat, 8" x 8", with logo

11210-000036

AED Location Sign

Flat, 10" x 10", with logo

11210-000037

AED Location Sign

Flat, 8" x 8", without logo

11210-000038

AED Location Sign

Flat, 10" x 10", without logo

11210-000039



AED Location Sign

T-mount, 8" x 8", with logo

11210-000031

AED Location Sign

T-mount, 10" x 10", with logo

11210-000032

AED Location Sign

T-mount, 8" x 8", without logo

11210-000040

AED Location Sign

T-mount, 10" x 10", without logo

11210-000041

LITERATURE AND TRAINING MATERIAL

LIFEPAK CR Plus/LIFEPAK EXPRESS AEDs Reference Manual CD

21300-006380

LIFEPAK CR Plus/LIFEPAK EXPRESS AEDs Owners Manual

26500-001999

LIFEPAK CR Plus/LIFEPAK EXPRESS AEDs Orientation DVD

26500-002274

AED CPR Reference Poster

Poster-sized easy reference AED and CPR instruction guide. Applicable for wall-mounting near an AED.

26500-000268

Quick Reference Instruction Card

Laminated easy reference AED and CPR instruction card. Fits inside the lid.

26500-001096

AED Challenge®, CD-ROM

An interactive refresher training program that simulates cardiac arrest scenarios using LIFEPAK automated external defibrillators. PC training tool provides practice for AED personnel who are just learning to use an AED, and helps experienced AED users keep their skills sharp. Consistent with the American Heart Association's 2005 Guidelines. Not a replacement for initial training. Designed for responders who have completed initial CPR/AED classroom training.

26996-000015

AED Challenge, Web-Enabled

An interactive refresher training program that simulates cardiac arrest scenarios using LIFEPAK automated external defibrillators. Web-based training tool provides practice for AED personnel who are just learning to use an AED, and helps experienced AED users keep their skills sharp. Consistent with the American Heart Association's 2005 Guidelines. Not a replacement for initial training. Designed for responders who have completed initial CPR/AED classroom training.

Includes online tracking, reporting, and administration tools. Annual license for individual users is required. List price is per user/per year. Quantity discounts are available.

Web version is currently available to Heart Safe Solutions customers only.

26996-000014



LIFEPAK CR Plus AED Training System

Includes 1 AED trainer, remote control and cable, simulated CHARGE-PAK charging unit, five (5) pair training electrodes, carrying case, removable shock button cover, operating instructions, and batteries for the trainer and remote control. Mimics semiautomatic and fully automatic operation.

11250-000073

Operating Instructions for the LIFEPAK CR Plus AED Training System.....26500-001156

Replacement Carrying Case11260-000014

Replacement Remote Control and Cable11250-000024



Adult AED QUIK-PAK Training Electrode Set

Consists of five (5) pair of reusable AED training electrodes, cable connector assembly, and reusable foil pouch. Reusable training electrodes good for 30 applications.

11250-000012

Replacement Training Electrodes

Five (5) pair of adhesive electrode portion only.

11250-000015

Infant/Child AED Training Electrode Set

Consists of five (5) pair of reusable AED training electrodes, cable connector assembly, and reusable foil pouch. Reusable training electrodes good for 30 applications.

11250-000045

Cable/connector Assembly and Reusable Foil Pouch for Infant/Child AED Training Electrodes

Cable, connector and foil pouch only. Adhesive portion of AED training electrodes sold separately.

11250-000043

Replacement Infant/Child AED Training Electrodes

Five (5) pair of reusable training electrodes. Adhesive electrode portion only. Cable/connector assembly and reusable foil pouch sold separately.

11250-000042



©2008 Physio-Control, Inc. LIFEPAK, LIFEPAK CR, and LIFEPAK EXPRESS are registered trademarks and CHARGE-PAK and QUIK-PAK are trademarks of Physio-Control, Inc. AMBU, Res-Cue Mask and Rescue-Key are registered trademarks of AMBU, Inc. AED Challenge is a registered trademark of Insight Instructional Media, LLC. Microsoft and Windows are registered trademarks of Microsoft Corporation. Specifications subject to change without notice. Not all products available outside the U.S. Please consult your LIFEPAK Defibrillator Operating Instructions or Declaration of Conformity for a list of approved accessories in your region.

MIN 3207698-001 / CAT 26500-002683