

## IP Ratings Dust and Water resistance chart

First number (Protection against solid objects)	Definition	Second number (Protection against liquids)	Definition
0	No protection	0	No protection
1	Protected against solids objects over 50mm (e.g. accidental touch by hands	1	Protected against vertically falling drops of water
2	Protected against solids objects over 12mm (e.g. fingers)	2	Protected against direct sprays up to 15° from the vertical
3	Protected against solids objects over 2.5mm (e.g. tools and wires)	3	Protected against direct sprays up to 60° from the vertical
4	Protected against solids objects over 1mm (e.g. tools, wires and small wires	4	Protected against sprays from all directions - limited ingress permitted
5	Protected against dust - limited ingress (no harmful deposit)	5	Protected against low pressure jets if water from all directions - limited ingress permitted
6	Totally protected against dust	6	Protected against strong jets of water e.g. for use on ship decks - limited ingress permitted
		7	Protected against the effects of temporary immersion between 15cm and 1m. Duration of test 30 minutes
		8	Protected against long periods of immersion under pressure

HeartSine samaritan PAD - IP56 rating, highest in the industry for Dust and Water resistance